



LEVO Dealer:	Date:
Dealer reference / Client name:	Anglikerstrasse 20 CH-5610 Wohlen Phone: +41 56 618 44 11 Fax: +41 56 618 44 10 e-mail: office@levo.ch
Order number:	



Standard configuration 60.001.6601

- Transport dimensions: width 45 cm, length 62 cm, height 43 cm (incl. footplate)
- M/E Wheelchair: Manual driving, seating and powered standing
- Seat plate and backrest upholstery
- Flip-up armrests
- Armrests height and angle adjustable
- Armrests round with built-in control switch
- Fold down backrest
- Standard colors: platinum, red, yellow, blue or black
- One-piece fixed footrest
- Pneumatic rear wheels and pump
- Batteries and battery charger
- Chest strap and knee support with wedge

Anthropometric Data



1 Hip width

Please fill in the correct box according to the measured hip width in a sitting position:

30 cm 35 cm

Seat width:

30 cm	Option
35 cm	Option

2 Upper leg length

There are two seat frame sizes available for the Kid & Junior: XS and S. The upper leg length determines which frame size and which armrest size to choose for the wheelchair.

Please fill in the measured upper leg length:

Upper leg length:
28 cm 45 cm

Example:

Upper leg length 28 cm
 - 4 cm
Seat upholstery depth 24.0 cm

Seat upholstery depth:

XS	24.0 cm	Option
	26.0 cm	Option
S	28.5 cm	Option
	30.5 cm	Option
	33.0 cm	Option
	35.0 cm	Option

3 Lower leg length

The lower leg length determines the footrest position and the seat height of the wheelchair. The seat height determines the wheel size and with that the total height of the wheelchair.

Please fill in the measured lower leg length:

Lower leg length:
19 cm 41 cm

Please check the desired seat height and footrest position, under consideration of the seat cushion thickness.

Example:

Lower leg length 36 cm
Seat cushion height - 2.5 cm
Footrest position 33.5 cm

Now you may choose between seat height 33, 40 or 44 cm.

Footrest position:

Kid 33 with adapter	19.2 cm	Option
	21.0 cm	Option
Junior w/o	22.8 cm	Option
	24.6 cm	Option
Junior 40 with adapter	26.4 cm	Option
	28.2 cm	Option
Junior 44 w/o adapter	30.0 cm	Option
	31.8 cm	Option
	33.6 cm	Option
	35.4 cm	Option
	35.2 cm	Option
	39.0 cm	Option
	40.8 cm	Option
	Junior 44 with adapter	

Footrest adapter for short legs (Junior)

60.003.6610

4 Elbow height

Standard armrest height is 16 cm

10 cm 22 cm

Please measure the anthropometric data as shown. If possible, the user should be in a sitting position and wearing shoes.

Wheels and position, floor to seat height in the front, measured without seat cushion

Rear wheel 16", castor 3"
 Rear wheel 20", castor 5"
 Rear wheel 22", castor 5"

Kid: 33 (30) cm	Option
Junior: 40 (37) cm	Option
Junior: 44 (41) cm	Option

Stand-up operation



Right	Standard
Left	Option

Control options



① Hand control	60.005.0040
② Dual switch for table tray control, <i>table tray not included</i>	60.007.0010

Frame color

Platinum	Standard
Red	Option
Yellow	Option
Blue	Option
Black	Option
Blue metallic	60.008.0020
Green metallic	60.008.0020
Red metallic	60.008.0020
Silver metallic	60.008.0020
Orange reflex	60.008.0030
Lemon reflex	60.008.0030
Silver reflex	60.008.0030
Aluminium polished	60.008.0050
Other RAL color	RAL: _____ 60.008.0010

Backrest




Standard sling backrest height 24 cm	Height: _____	Standard
Backrest extension 15 cm, <i>for backrest height 30 cm the extension must be added</i>		32.066.0410
① V-Trak backrest, height 28 cm, <i>only for seat width 30 cm</i>		60.009.6830
V-Trak backrest, height 30 cm		60.009.6610
② V-Trak Axxis reduced height backrest, height 31 cm		60.009.0070
V-Trak Axxis low backrest, height 40 cm		60.009.0050
V-Trak backrest extension, 13 cm		32.066.0480
<i>without backrest, be aware net price is lower after discount</i>		60.109.6810


Backrest angle

85°	Option
90°	Standard
95°	Option

Seat cushion with memory foam depending on the seat depth and seat width




Seat width ►	30	35
Seat depth ▼		
XS	30.100.3760	30.100.0610
S	30.100.0660	30.100.3740




Seat cushion cover

Seat width ►	30	35
Seat depth ▼		
XS	30.100.3761	30.100.0611
S	30.100.0661	30.100.3741




Urin proof seat cushion cover

Seat width ►	30	35
Seat depth ▼		
XS	on request	on request
S	30.100.0663	30.100.3743



Headrest



① V-Trak headrest standard, <i>only in combination with V-Trak Axxis backrest</i>	02.067.116
② V-Trak headrest winged, <i>only in combination with V-Trak Axxis backrest</i>	02.067.117

Lateral supports and hardware



V-Trak lateral support, pad and cover

Medium	Right		32.068.0030
Medium	Left		32.068.0040

V-Trak lateral support, pad and cover

Large	Right		32.068.0010
Large	Left		32.068.0020

Knee support



① Rectangular knee wedge

50x100		Option
60x100		Option
50x150		Option

② Conical knee wedge, *gives more support, in particular with spasticity*

50-100x100		60.011.0020
90-140x100		60.011.0020
130-180x100		60.011.0020

Foot plate

Foot plate angle

70°		Option
80°		Option
90°		Standard
100°		Option
110°		Option

Non-slip footrest surface

60.003.0040

Wheel options



① Spinegy wheels, *only for 22" wheels*

Red		60.002.0015
Black		60.002.0015
White		60.002.0015


② Maxgrepp original hand ring, with rubber cover, *only for 22" wheels*

60.002.0060

③ Maxgrepp Ergo Para CP hand ring, *only for 22" wheels*

60.002.0050

Accessories



① Table tray

32.070.0115

② Folding down stroller bar, *width adjustable*

02.069.0030

③ Anti-tip for Kid

Right		32.066.0420
Left		32.066.0420

③ Anti-tip for Junior

Right		32.066.0430
Left		32.066.0430

④ Spoke guard transparent

Paar	20"		02.066.1013
Paar	22"		02.066.1014

⑤ Backpack

Red		30.100.0841
Yellow		30.100.0842
Blue		30.100.0843
Black		30.100.0844

⑥ Foot guide

Right		32.069.95
Left		32.069.96

⑦ Brake lever extension (6 cm)

Right		30.100.5030
Left		30.100.5030

⑧ Knee support height adapter, *for users with long lower leg lengths, knee holder gets closer to the joint*

Pair 32.069.85

⑨ Control switch extension

Vertical		32.066.0450
Horizontal		32.070.0156

⑩ Quick adjustment for seat depth / footrest (Set)
Quick adjustment for seat depth (Set)

60.004.0020
60.004.0030

