

LEVO Dealer:

Date:

Dealer reference / Client name:

Anglikerstrasse 20  
CH-5610 Wohlen  
Phone: +41 56 618 44 11  
Fax: +41 56 618 44 10  
e-mail: office@levo.ch

Order number:


Standard configuration

60.001.7201

- Power Base: width 63 cm; length 105 cm (incl. footplate)
- E/E Wheelchair with 4 active wheels and "Low shearing system"
- Seat support plate
- Seat height 48 cm
- V-Trak Axxis low back, 40 cm/16" high
- Quick release for reclining back
- Flip-up armrests
- VR2 control system from PG Drives Technology
- 48 Ah Batteries and battery charger
- Standard colors: red, yellow, blue or black
- Flip-up footrest side- and backwards with non slip surface
- Chest strap and knee support pro
- 8" single rear wheel castors
- Vehicle tie-down locations



### Anthropometric Data



Please measure the anthropometric data as shown. If possible, the user should be in a sitting position and wearing shoes.

<b>1 Hip width</b>															
Please check the correct box according to the measured hip width in a sitting position:															
32 cm	<input type="checkbox"/> 52 cm														
<table border="1"> <tr> <td colspan="2"><b>Seat width:</b></td> </tr> <tr> <td>32 cm</td> <td><input type="checkbox"/> Option</td> </tr> <tr> <td>36 cm</td> <td><input type="checkbox"/> Option</td> </tr> <tr> <td>40 cm</td> <td><input type="checkbox"/> Option</td> </tr> <tr> <td>44 cm</td> <td><input type="checkbox"/> Option</td> </tr> <tr> <td>48 cm</td> <td><input type="checkbox"/> Option</td> </tr> <tr> <td>52 cm</td> <td><input type="checkbox"/> Option</td> </tr> </table>		<b>Seat width:</b>		32 cm	<input type="checkbox"/> Option	36 cm	<input type="checkbox"/> Option	40 cm	<input type="checkbox"/> Option	44 cm	<input type="checkbox"/> Option	48 cm	<input type="checkbox"/> Option	52 cm	<input type="checkbox"/> Option
<b>Seat width:</b>															
32 cm	<input type="checkbox"/> Option														
36 cm	<input type="checkbox"/> Option														
40 cm	<input type="checkbox"/> Option														
44 cm	<input type="checkbox"/> Option														
48 cm	<input type="checkbox"/> Option														
52 cm	<input type="checkbox"/> Option														
<b>2 Upper leg length</b>															
Please fill in the measured upper leg length:															
38 cm	<input type="checkbox"/> 66 cm														
<b>3 Lower leg length</b>															
Please fill in the measured lower leg length:															
24 cm	<input type="checkbox"/> 44 cm														
If the lower leg length is between 45 and 50 cm, the C3 is available with a 6 cm higher seat height. The seat height is not adjustable. Please check the box accordingly.															
	<input type="checkbox"/> Option														
<b>4 Elbow height</b>															
Standard armrest height is 22 cm															
17 cm	<input type="checkbox"/> 26 cm														
<b>User weight</b>															
The weight of the user determines the force of the gas springs.															
The maximum user weight is 140 kilograms. Please fill in the box.															
0 - 60 kg	1200 N														
60 - 100 kg	1700 N <input type="checkbox"/> kg														
100 - 140 kg	2100 N														

### Controls



<b>1</b> VR2 Control system: - Standard control system - 90 A Power module - Max. 2 actuators	- Controls lights - Programmable	<input type="checkbox"/>	Standard
<b>2</b> R-Net control system "Standard": - Color LCD Joystickmodule - 120 A Power module - Max. 6 actuators, multiple axes	- Controls lights - Fully programmable, multiple profiles - Expandable	<input type="checkbox"/>	60.072.0010
<b>3</b> R-Net control system "XL", as R-net control system "Standard" with the following extras : - 3.5" VGA color display - Paddle switches for on/off, profile/mode and speed	- Metal case - Loudspeaker instead of a buzzer for a louder horn	<input type="checkbox"/>	60.072.0020
<b>3</b> R-Net control system "XL", with integrated bluetooth		<input type="checkbox"/>	60.072.0030

**Joystick position**



- ① Swing-away, offers more space for transfers and allows to drive closer to tables
- ② Joystick protection guard, shock protection for joystick, protects the users hands

<b>Right</b>	Standard
<b>Left</b>	Option
	60.005.7210
	60.012.6860

**Attendant Controls**



VR2 Attendant control  
 R-Net Attendant control

For special controls, see our special controls order form.

	32.072.0140
	32.068.1850

**Speed**

6 km/h  
 8 km/h  
 10 km/h

	Option
	Option
	Standard

**Stability\***



Gyro module, only available for R-net control system

	32.072.0400
--	-------------

- Perfect control on cambers
- Straight line driving over obstacles
- Eliminates spin-out and fish-tailing

\* This option is only available in Europe

**Batteries**

48 Ah  
 73 Ah

Without batteries, be aware net price is lower after discount  
 Without charger, be aware net price is lower after discount

	Standard
	60.006.7210
	60.106.6810
	60.111.6710

**Frame color**

Red  
 Yellow  
 Blue  
 Black  
 Blue metallic  
 Green metallic  
 Red metallic  
 Silver metallic  
 Orange reflex  
 Lemon reflex  
 Silver reflex  
 Traffic white  
 Other RAL color

RAL: \_\_\_\_\_


	Standard
	Option
	Option
	Option
	60.008.0020
	60.008.0020
	60.008.0020
	60.008.0020
	60.008.0030
	60.008.0030
	60.008.0030
	60.008.0030
	60.008.0010
	60.008.0010

**Tube color**

Standard grey metallic  
 Matt black

	Standard
	Option

**Wheels**



- ① Block profile drive wheels
  - ② V-profile drive wheels
- Puncture proof tires on drive wheels  
 Air tires on drive wheels
- Black tires
- ③ 7" twin rear castors

	Standard
	Option
	Standard
	Option
	On request
	60.002.0140

**Lights**



Lights and indicators for outdoor use

	60.012.7220
--	-------------

**Seating**

**Electrical adjustments**




Tilt in space 0 - 35°

Reclining backrest, only in combination with R-Net

	60.012.7210
	60.009.0150

**Backrest, The V-Trak Axxis low back is delivered with the chair as a standard feature**



- ① V-Trak Axxis low backrest, height 40 cm/16"
- ② V-Trak Axxis high backrest, height 61 cm/24"
- V-Trak Axxis reduced height backrest, height 31 cm/12"

	Standard
	60.009.0080
	60.009.0070

- ③ V-Trak Axxis Deep contour

Size	M	L	XL	XXL
Type	inside 26 - 37 cm	inside 30 - 40 cm	inside 33 - 46 cm	inside 35 - 52 cm
a High				
b Reduced height				

	60.009.0110
	60.009.0130

- ④ Third party backrest conversion bracket, *please choose the correct size*

32 - 36	32.072.0351
40 - 44	32.072.0352
44 - 52	32.072.0353

Without backrest shell, *be aware net price is lower after discount*

	60.109.7210
--	-------------

**Seat cushion, with memory foam depending on the seat depth and seat width incl. cover**



Seat width	32	36	40	44	48	52
35, 37, 39 cm	30.103.5531	30.103.5532	30.903.5535	30.903.5531	30.903.5538	on request
41, 43, 45 cm	30.903.5532	30.103.5533	30.103.5534	30.103.5539	30.903.5536	30.903.5537
47, 49, 51 cm	30.903.5533	30.103.5535	30.103.5536	30.103.5541	30.103.5542	30.103.5549
53, 55, 57 cm	on request	30.903.5534	30.103.5537	30.103.5543	30.103.5544	30.103.5545
59, 61, 63 cm	on request	on request	30.103.5538	30.103.5546	30.103.5547	30.103.5548

**Seat cushion cover**

Seat width	32	36	40	44	48	52
35, 37, 39 cm	30.103.6211	30.103.6212	30.903.6215	30.103.6220	on request	on request
41, 43, 45 cm	30.903.6212	30.103.6213	30.103.6214	30.103.6219	30.103.6235	30.103.6230
47, 49, 51 cm	30.903.6213	30.103.6215	30.103.6216	30.103.6221	30.103.6222	30.103.6229
53, 55, 57 cm	on request	30.903.6214	30.103.6217	30.103.6223	30.103.6224	30.103.6225
59, 61, 63 cm	on request	on request	30.103.6218	30.103.6226	30.103.6227	30.103.6228

Seat cushion cover urin proof for all seat depths and widths

	30.903.6230
--	-------------

**Headrest**



- ① V-Trak headrest standard, *only in combination with V-Trak Axxis backrest*
- ② V-Trak headrest winged, *only in combination with V-Trak Axxis backrest*
- Track extension for V-Trak head support

	02.067.116
	02.067.117
	02.009.0028

**Knee support**



- ① Knee support Pro
- ② Knee support Integral
- Extension for Integral +100 mm

Size S (10 - 12 cm)	Standard
Size M (14 - 16 cm)	Standard
Size S (10 - 12 cm)	60.011.0040
Size M (14 - 16 cm)	60.011.0040
	Option

**Accessories**



- ① Skirt guard
- ② Upper leg support
- Holder to skirt guard and upper leg support

Right	32.072.0021
Left	32.072.0022
Right	32.072.0031
Left	32.072.0032

32	36	40	44	48	52
32.072.0151	32.072.0152	32.072.0153	32.072.0154	32.072.0155	32.072.0156

**Accessories**



③ Chest role

32	36	40	44	48	52
32.072.0041	32.072.0042	32.072.0043	32.072.0044	32.072.0045	32.072.0046

④ Foot guide

<b>Small</b>	<b>Right</b>		32.069.95
<b>Small</b>	<b>Left</b>		32.069.96
<b>Large</b>	<b>Right</b>		32.070.0040
<b>Large</b>	<b>Left</b>		32.070.0041

⑤ Push handle bar

	32.072.0090
--	-------------

⑥ Rearview mirror

<b>Right</b>	32.072.0050
<b>Left</b>	32.072.0050

⑦ Swing-away table tray

<b>32 - 40</b>	32.072.0011
<b>44 - 52</b>	32.072.0012

Mounted on the opposite side of the joystick position

<b>Right</b>	
<b>Left</b>	

Quick release kit for seat depth and footrest position adjustment

	60.004.7200
--	-------------

⑧ Joystick top

<b>T-bar</b>	20.067.22
<b>Stick</b>	20.067.23
<b>Softball</b>	20.067.24
<b>Mushroom</b>	20.067.25
<b>Chin</b>	20.067.26
<b>Ball</b>	20.067.27

Remarks