

LEVO Dealer:

Date:

Dealer reference / Client name:

Anglikerstrasse 20  
CH-5610 Wohlen  
Phone: +41 56 618 44 11  
Fax: +41 56 618 44 10  
e-mail: office@levo.ch

Order number:


Standard configuration

60.001.6801

- Power Base: width 66 cm, length 110 cm (incl. footplate)
- E/E Wheelchair: powered driving, sitting and standing
- Seat plate and backrest upholstery
- Flip-up armrests
- Height adjustable armrests
- Backrest length adjustment and armrests height adjustment
- Manual seat height adjustment
- R-Net control system with color display
- 48 Ah batteries and charger
- Standard colors: white, red, blue or black
- Flip-up footrest
- Chest strap and knee support with wedge



### Anthropometric Data



#### 1 Hip width

Please check the correct box according to the measured hip width in a sitting position:

33 cm  50 cm

Seat width:

38 cm	<input type="checkbox"/>	Option
44 cm	<input type="checkbox"/>	Option
50 cm	<input type="checkbox"/>	Option

#### 2 Upper leg length

Please fill in the measured upper leg length, seat depth is 3 cm shorter.

35 cm  59 cm

Seat depth:

S-M Frame	36.0 cm	<input type="checkbox"/>	Option
	38.2 cm	<input type="checkbox"/>	Option
	40.4 cm	<input type="checkbox"/>	Option
	42.6 cm	<input type="checkbox"/>	Option
	44.8 cm	<input type="checkbox"/>	Option
L-XL Frame	47.0 cm	<input type="checkbox"/>	Option
	49.2 cm	<input type="checkbox"/>	Option
	51.4 cm	<input type="checkbox"/>	Option
	53.6 cm	<input type="checkbox"/>	Option
	55.8 cm	<input type="checkbox"/>	Option
	58.0 cm	<input type="checkbox"/>	Option

#### 3 Lower leg length

Please fill in the measured lower leg length:

37 cm  64 cm

The seat height is dependent on the lower leg length. With a longer lower leg length the seat is higher.

Lower leg length in cm

Seat height in cm

37,0-44,2	46,0-47,8	49,6-51,4	53,2-55,0	56,8-58,6	60,4-62,2	64,0
56	58.5	61	63.5	66	68.5	71

#### 4 Elbow height

Standard armrest height is 20 cm

17 cm  26 cm

Backrest height, for standard backrest

Please measure the distance between seat plate and lowest shoulder blade point:

cm

User weight

The weight of the user determines the force of the support wheel gas springs.  
The maximum user weight is 120 kilograms. Please fill in the box.

kg

Please measure the anthropometric data as shown. If possible, the user should be in a sitting position and wearing shoes.

**Controls**



①	R-Net control system "Standard": - Color LCD Joystickmodule - 120 A Power module - Max. 6 actuators, multiple axes - Controls lights - Fully programmable, multiple profiles - Expandable		Standard
②	R-Net control system "XL", as R-net control system "Standard" with the following extras : - 3.5" VGA color display - Paddle switches for on/off, profile/mode and speed - Metal case - Loudspeaker instead of a buzzer for a louder horn		60.075.0010
②	R-Net control system "XL" with integrated bluetooth		60.075.0030

**Joystick position**



		<b>Right</b>	Standard
		<b>Left</b>	Option
①	Swing-away, offers more space for transfers and allows to drive closer to tables		60.005.6810
②	Joystick protection guard, shock protection for joystick, protects the users hands		60.012.6860

**Attendant Controls**



R-Net Attendant control		32.068.1850
Connector Block for R-Net, needed for special controls		32.068.1880
For special controls, see our special controls order form.		

**Speed**

6 km/h		Option
8 km/h		Option
10 km/h		Standard

**Stability\***



Gyro module		32.068.1000
- Perfect control on cambers - Straight line driving over obstacles - Eliminates spin-out and fish-tailing		
* This option is only available in Europe		


**Batteries**

48 Ah		Standard
Without batteries, be aware net price is lower after discount		60.106.6810
Without charger, be aware net price is lower after discount		60.111.6710

**Frame color**

White		Standard
Red		Option
Blue		Option
Black		Option
Blue metallic		60.008.0020
Green metallic		60.008.0020
Red metallic		60.008.0020
Silver metallic		60.008.0020
Orange reflex		60.008.0030
Lemon reflex		60.008.0030
Silver reflex		60.008.0030
Other RAL color	<b>RAL:</b> _____	60.008.0010

**Wheels**



①	Block profile drive wheels		Standard
②	V-profile drive wheels		Option
③	Curb Climber		60.002.0100

**Lights**



Lights and indicators for outdoor use		60.012.6810
---------------------------------------	--	-------------


**Seating**



**Electrical adjustments**

- Powered seat lift 150 mm 60.012.6820
  - ① Powered leg rest, *makes an angle adjustment of 90-180° possible, independent from the backrest* 60.012.6870
  - ② Combination of powered reclining back and powered leg rest  
*Enables a complete lying position and the function lying to standing is integrated. Enables a lying transfer.* 60.012.6830
- Recommendation: order with additional V-Trak arms and head rest.*

**Backrests**



- Standard backrest upholstery Standard
  - V-Trak Axxis reduced height backrest, height 31 cm 60.009.0070
  - V-Trak Axxis low backrest, height 40 cm 60.009.0050
  - V-Trak Axxis high backrest, height 61 cm 60.009.0060
  - ③ V-Trak Axxis Deep contour
- | Size             | M                 | L                 | XL                | XXL               |
|------------------|-------------------|-------------------|-------------------|-------------------|
| Type             | inside 26 - 37 cm | inside 30 - 40 cm | inside 33 - 46 cm | inside 35 - 52 cm |
| a High           |                   |                   |                   |                   |
| b Reduced height |                   |                   |                   |                   |
- 60.009.0120  
60.009.0140
- Without backrest, *be aware net price is lower after discount* 60.109.6810

**Backrest options**



- ① Backrest extension 15 cm, *for standard backrest* **38**  
**44**  
**50** 32.062.3815  
32.062.4415  
32.062.5015
  - ② Backrest extension with push handle bar and upholstery, *for standard backrest*
- | Seat width | 38          | 44          | 50          |
|------------|-------------|-------------|-------------|
| Seat depth |             |             |             |
| 15 cm      | 32.068.0761 | 32.068.0771 | 32.068.0781 |
| 25 cm      | 32.068.0762 | 32.068.0772 | 32.068.0782 |
| 35 cm      | 32.068.0763 | 32.068.0773 | 32.068.0783 |
- V-Trak backrest extension, 20 cm/8" **38 - 50** 32.069.99
- Additional V-Trak arms, *recommended when selecting reclining back* 60.009.6810

**Seat cushion with memory foam depending on the seat depth and seat width**



Seat width	38	44	50
Seat depth			
S	30.102.9273	30.102.9274	30.102.9275
M	30.102.9276	30.102.9277	30.102.9278
L	30.102.9281	30.102.9282	30.102.9283
XL	30.102.9284	30.102.9285	30.102.9286

**Additional seat cushion cover**

Seat width	Standard			Urin Proof		
	38	44	50	38	44	50
Seat depth						
S	30.102.9461	30.102.9462	30.102.9463	30.102.9481	30.102.9482	30.102.9483
M	30.102.9464	30.102.9465	30.102.9466	30.102.9484	30.102.9485	30.102.9486
L	30.102.9467	30.102.9468	30.102.9469	30.102.9487	30.102.9488	30.102.9489
XL	30.102.9471	30.102.9472	30.102.9473	30.102.9474	30.102.9475	30.102.9476

**Headrest**



- ① V-Trak headrest standard, *only in combination with V-Trak and Axxis backrest* 02.067.116
  - ② V-Trak headrest winged, *only in combination with V-Trak and Axxis backrest* 02.067.117
- Track extension for V-Trak head support, for V-Trak 60 and 80 cm 02.009.0028

**Lateral supports and hardware**



- V-Trak lateral support, pad and cover **Medium Right**  
**Medium Left** 32.068.0030  
32.068.0040
- V-Trak lateral support, pad and cover **Large Right**  
**Large Left** 32.068.0010  
32.068.0020

## Knee support



Standard knee support

	Standard
--	----------

- ① Rectangular knee wedge

50 x 100	Option
60 x 100	Option
80 x 100	Option
50 x 150	Option
80 x 150	Option



- ② Conical knee wedge, gives more support, in particular with spasticity

50-100 x 100	60.011.0020
90-140 x 100	60.011.0020
130-180 x 100	60.011.0020



- ③ Knee support "Pro"

Foam size S (10 - 12 cm)	60.011.0030
Foam size M (14 - 16 cm)	60.011.0030

LEVO®

## Foot plate

Foot plate depth

27 cm	Standard
20 cm	60.003.0090

Foot plate width, 30 cm only possible with lower leg lengths between 23 and 27 cm

22 cm	Standard
30 cm	60.003.0110

Foot plate angle, foot plate is angle adjustable in all three positions

80°	Option
90°	Standard
100°	Option



- ① Foot plate wedge, better positioning with spasticity or foot deformation

50-100 x 80	30.101.8171
80-130 x 80	30.101.8172



- ② Foot protection frame, protects the users feet

	60.003.0100
--	-------------

Non-slip footrest surface


	02.001.0702
--	-------------

## Accessories



- ① Skirt guard

S - M	32.068.0830
L - XL	32.068.0840



- ② Upper leg support

Right	32.068.1751
Left	32.068.1752



- ③ Chest role

38	44	50
32.068.1761	32.068.1762	32.068.1763



- ④ Foot guide

Small	Right	32.069.95
Small	Left	32.069.96
Large	Right	32.070.0040
Large	Left	32.070.0041



- ⑤ Push handle bar, can be used with Standard and V-Trak backrest

38	44	50
32.068.0791	32.068.0792	32.068.0793



- ⑥ Rearview mirror

Right	32.068.0120
Left	32.068.0120



- ⑦ Swing-away table tray, can only be mounted on the opposite side of the joystick position

38	44	50
32.068.1061	32.068.1062	32.068.1062

Right	
Left	

## Accessories



- ⑧ Fixed table tray

38	44	50
32.062.0038	32.062.0044	32.062.0050

Vehicle tie-down kit

	32.068.1612
--	-------------

- ⑨ Support for QLK docking system, without assembly cost  
Assembly cost for EZ lock brackets

	32.068.2050
	33.000.0000

Cover against clamping danger

	Option
--	--------

**Accessoires**



⑩ Joystick top

<b>T-bar</b>		20.067.22
<b>Stick</b>		20.067.23
<b>Softball</b>		20.067.24
<b>Mushroom</b>		20.067.25
<b>Chin</b>		20.067.26
<b>Ball</b>		20.067.27

Remarks